SPENCER®

Rigid splints with flexible internal structure Stecchobende rigide con anima flessibile

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User's Manual

This appliance conforms with the Directive 93/42/CEE "Medical Devices".
Guarantee of Quality system for the production and the final control of the products certified by the notifying body TÜV Product Service GmbH.

1 GENERAL INFORMATION

1.1 Aim and contents

The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.

1.2 Conservation of the instruction manual

The instruction and maintenance manual must be kept with the product, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.

1.3 Symbols used

General or specific warning

See instructions for use

Lot number

Product code

Product compliant with specifications of the Directive 93/42/CEE

1.4 Servicing request

For any information regarding the use, maintenance and installation, please contact the Spencer Customer Care Service on tel. 0039 0521 541111, fax 0039 0521 541222, e-mail info@spencer.it or write to Spencer Italia S.r.l. - Strada Cavi, 7 - 43044 Collecchio (Parma) - ITALY.

1.5 Demolition

Follow the current regulations.

1.6 Labelling

The serial number as indicated below can be found on each appliance and must not be removed or covered. In order to facilitate assistance please indicate or communicate the lot number (LOT) on the label.

2 WARNINGS

2.1 General warnings

* Before carrying out any kind of operation on the appliance, the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.

* In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.

* Regularly check the appliance.

* In the case of any abnormalities or damage to the appliance, which could jeopardize the functioning and the safety, the appliance must be immediately removed from service.

* Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or resuer.

* The appliance must not in any way be tampered with. In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself.

* Those who modify or have modified, prepared or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.

* Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.

* Handle with care.

2.2 Specific warnings

* The product must be used by trained personnel only.

* If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.

* When the device is being used, the assistance of qualified staff must be guaranteed.

* The splint must be applied by at least two trained operators, with good practical sense and common sense.

* Before modelling the splint, evaluate the application in relation to the visibility of wounds or fractures.

* The device should not be exposed nor get in contact with heat sources or flammable agents.

* Use the device only as described in this manual.

* Always check the integrity of all the parts of the device before use.

2.3 Contraindications and side effects

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

3 PRODUCT DESCRIPTION

3.1 Intended use

Blue/Splint/Blue/Splint Pro are flexible splints designed to immobilise dislocated or sprained limbs. They can be adapted to the various positions indicated for the pathology and are quick to apply. They are not bulky once applied and they do not increase the patient's inertial mass. Their flexibility is due to the possibility of shaping the metal core inside the Blue/Splint/Blue/Splint Pro in order to adapt it to the shape of the limb to immobilize in order to position the limb in a reliable manner for stability during transport. Each splint can therefore be adapted to the fracture type and the rescuer does not have to bother to reduce the dimensions. This permits transport of injured limbs in angled pathological position which will guarantee the correct limitation of movement in all directions. The closures in Velcro® (for Blue/Splint) or with the rapid adjustment and closure system (Blue/Splint Pro) are quick to apply and are easily adapted to the various patient sizes. Use of the various sizes offers the possibility to immobilise multiple fractures. Stored in their special transport bag, or applied to the patient, the Blue/Splint/Blue/Splint Pro are characterised by their particularly reduced volume which makes them ideal for any type of operation in any condition. The patient cannot intervene on the device in any way.

3.2 Main components (fig. 8)

A Nyprene structure (the outside part is blue, the inside is black)

B Velcro® closures (the orange webbing highlights the top part, the blue webbing is used for the other points)

C Quick-release and regulation closures (nylon webbing with integrated regulator)

D Flexible inside structure made of a tempered ergal profile

E Ripstop carrying case

3.3 Models

JM80003A Blue Splint leg

JM80004A Blue Splint arm

JM80005A Blue Splint forearm

JM80006A Blue Splint wrist

JM80007A Blue Splint elbow/ankle

JM80002A Blue Splint Kit 5 sizes with bag

JM80030A Blue Splint Pro leg

JM80031A Blue Splint Pro arm

JM80032A Blue Splint Pro forearm

JM80033A Blue Splint Pro wrist

JM80034A Blue Splint Pro elbow/ankle

JM80035A Blue Splint Pro Kit 5 sizes with bag

3.4 Technical data

Dimensions Blue/Splint/Blue/Splint Pro leg: 578 x 420 mm

Dimensions Blue/Splint/Blue/Splint Pro arm: 582 x 290 mm

Dimensions Blue/Splint/Blue/Splint Pro forearm: 387 x 300 mm

Dimensions Blue/Splint/Blue/Splint Pro wrist: 302 x 255 mm

Dimensions Blue/Splint/Blue/Splint Pro elbow/ankle: 540 x 305 mm

Dimensions Blue/Splint/Blue/Splint Pro Kit: 640 x 100 x h310 mm

Weight Blue/Splint/Blue/Splint Pro Kit: 680 x 80 x h350 mm

Weight Blue/Splint/Blue/Splint Pro Kit: 1.5 kg

Functioning temperature: from -40 to +60 °C

Storage temperature: from -40 to +86 °C

Relative humidity: from 30 to 75%

Blue/Splint/Blue/Splint Pro are manufactured in a material which permits excellent thermal isolation of the immobilised limb and they permit the application of localised pressure on the fixing points. The use of Velcro® (or of the rapid adjustment and attachment system) with the elastic properties of the Nyprene material permit complete control of compression aided by the elastic properties of the material. All materials used are long lasting and 100% impermeable.

4 OPERATING INSTRUCTIONS

4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packed ensuring also that there are no risks of shocks, bumps or falls during the transport itself. Keep the original packaging for use in case of any further transport. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The appliance must be stored in a dry place free from humidity.

4.3 Functioning

See back.

4.4 Troubleshooting

Problem

1. Closure is insufficient

2. The metal core inside has lost its functional characteristics

Cause

1. Closing system is worn out

2. Wear and tear

Remedy

1-2. Put immediately the device out of service and contact the service centre

5 MAINTENANCE AND CLEANING

5.1 Cleaning

Can be carried out using alcohol or other neutral disinfectants. If the metal core is removed, they can be machine washed at the maximum temperature of 40 °C. The straps of the Blue/Splint Pro model must be closed before machine washing.

5.2 Maintenance

The person responsible for every day maintenance can only substitute the spare parts indicated on paragraph 6.2 "Spare Parts". All other substitutions or repairs can be carried out only by the manufacturer or by a centre authorised by the manufacturer.

The device does not require programmed servicing.

5.2.2 Special servicing
Only the manufacturer or centres with written authorisation are authorised to complete any special servicing operations. For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device. The device, if used as indicated in the following instruction manual, has an average life span of 5 years. The life span can be expanded only following a general revision of the product that must be carried out by the manufacturer or by a centre authorised by the manufacturer.

6 ACCESSORIES AND SPARE PARTS

6.1 Accessories

JM80008A Ripstop carrying case

6.2 Spare parts

JM80014A Metal profile for arm splint

JM80015A Metal profile for forearm splint

JM80016A Metal profile for wrist splint

JM80017A Metal profile for elbow/ankle splint

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